




Liquid Handling Consumables

Fisherbrand Elite Pipette and SureOne Pipette Tips



Fisherbrand™ Elite™ Pipettes and Fisherbrand™ SureOne™ Pipette Tips work together to form a complete, high-quality pipetting system.

Elite Pipettes

Achieve the best results in your lab with exceptional ergonomics while excelling at accuracy and precision. This comprehensive product line offers a wide range of pipettes for multiple applications, providing an optimal fit with SureOne Pipette Tips.

Elite Pipettes have an ergonomic design with index finger and hand support for maximum comfort. The low plunger force and larger soft-touch tip ejection helps you work longer without fatigue.

The entire pipette can be steam autoclaved to ensure sterility. Its double-action pipetting button has a rotating upper disc that prevents accidental volume changes during use, and a large display lets you see your selected volume more easily. These pipettes are available in a range of volumes that offer accuracy and precision for all your applications.



- ① Comfortable, lightweight handle with finger rest
- ② Index finger support
- ③ Hand support
- ④ Optimized fit for SureOne Tips
- ⑤ Soft-touch tip ejector
- ⑥ Extremely large display
- ⑦ Large tip ejector
- ⑧ Double-action pipetting button



Variable Volume Single-Channel Pipettes

IBS Code	Range	Increment	Accuracy	Precision	SureOne Tip for Calibration
11815762	0.2 to 2µL	0.002µL	2.5 to 12.0%	2.0 to 10%	11953416
11825762	0.5 to 5µL	0.01µL	1.5 to 6.0%	1.0 to 5.0%	11953416
11835762	1 to 10µL	0.02µL	1.0 to 2.5%	0.5 to 2.0%	11953416
11845762	2 to 20µL	0.02µL	1.0 to 3.0%	0.4 to 1.5%	10678325
11855762	5 to 50µL	0.1µL	0.6 to 3.0%	0.3 to 2.5%	10678325
11865762	10 to 100µL	0.2µL	0.8 to 3.0%	0.2 to 1.0%	10678325
11875762	20 to 200µL	0.2µL	0.6 to 1.8%	0.2 to 0.7%	10678325
11885762	100 to 1000µL	1µL	0.6 to 1.0%	0.2 to 0.6%	11568442
11895762	0.5 to 5mL	0.01mL	0.5 to 2.0%	0.2 to 0.8%	11917754
11805772	1 to 10mL	0.02mL	0.5 to 2.0%	0.2 to 0.8%	11957754

IBS Code	Description
15268638	Fisherbrand Elite Kit, Includes Stand and Four Pipettes: 0.2 to 2µL, 2 to 20µL, 20 to 200µL, 100 to 1000µL

Variable Volume Multichannel Pipettes

IBS Code	Channels	Range	Increment	Accuracy	Precision	SureOne Tip for Calibration
11815772	8	1 to 10µL	0.02µL	2.4 to 12%	1.6 to 8.0%	11953416
11835772	8	5 to 50µL	0.1µL	1.5 to 5.0%	0.7 to 2.0%	10678325
11825772	8	10 to 100µL	0.2µL	1.3 to 5.0%	0.5 to 2.0%	10678325
11845772	8	30 to 300µL	1.0µL	1.0 to 4.6%	0.3 to 1.5%	11903456
11855772	12	1 to 10µL	0.02µL	2.4 to 12%	1.6 to 8.0%	11953416
11875772	12	5 to 50µL	0.1µL	1.5 to 5.0%	0.7 to 2.0%	10678325
11865772	12	10 to 100µL	0.2µL	1.3 to 5.0%	0.5 to 2.0%	10678325
11885772	12	30 to 300µL	1.0µL	1.0 to 4.6%	0.3 to 1.5%	11903456

IBS Code	Description
11895772	Fisherbrand Elite Stand, Holds Six Pipettes

SureOne Pipette Tips

The line of universal fit pipette tips is available in many packaging configurations and volumes from 5µL to 10mL.

- Sterile tips are certified free of RNase, DNase, DNA, endotoxin, bioburden, PCR Inhibitor, and ATP
- Non-sterile tips are certified free of RNase, DNase, DNA, pyrogen, and endotoxin
- Non-sterile tips can be autoclaved at 120°C, 15 PSI or 1.034 bar for 15 minutes; do not autoclave sterile products
- Optimal fit with Fisherbrand™ Pipettes as well as universal fit on many other pipettes
- Standard and low-retention tips offered in filtered and non-filtered options
- Specialty tips like wide-orifice genomic tips, gel loading, and extended-length tips are available

SureOne Sterile Filtered Pipette Tips

SureOne Filter Tips use an inert, hydrophobic small micron filter that offers filter characteristics of aerosol resistance with optimized porosity, back pressure, and compression for maximum air flow and sample uptake. SureOne provides sample integrity through sample protection.

Volume	Style	Graduated	Tip Length	Color	Low Retention	Packaging	Part code	Fit
0.1 to 5µL	Micropoint Nano	At 2.5 µL	31 mm	Clear	No	Racked	11917734	Universal
					Yes	Reload *	11997714	Universal
0.1 to 10µL	Micropoint Micro	At 2µL	31 mm	Clear	No	Racked	11903466	Universal
					Yes	Racked	16490865	Universal
0.1 to 10µL	Micropoint Extended Length	At 2, 5, 10, and 20µL	46 mm	Clear	No	Racked	11913466	Universal
					Yes	Racked	16470865	Universal
0.1 to 20µL	Micropoint Extended Length	At 2, 5, 10, and 20µL	46 mm	Clear	No	Racked	11917724	Universal
					Yes	Racked	11923466	Eppendorf Style Pipettes
2 to 20µL	Beveled	No	51 mm	Clear	No	Racked	11977714	Universal
					Yes	Racked	16480865	Universal
2 to 20µL	Beveled	No	51 mm	Clear	No	Racked	11927724	Universal
					Yes	Reload *	11927724	Universal
2 to 20µL	Beveled	No	51 mm	Clear	No	Racked	11943466	Universal
					Yes	Racked	16420875	Universal
2 to 20µL	Beveled	No	51 mm	Clear	No	Reload **	11937724	Universal
					Yes	Racked	16430875	Universal
10 to 100µL	Beveled	No	51 mm	Clear	No	Racked	11953466	Universal
					Yes	Racked	16410875	Universal
10 to 100µL	Beveled	No	51 mm	Clear	No	Reload **	11947724	Universal
					Yes	Racked	11947724	Universal

*Use empty rack for 10 to 20µL: **11973456**

Use empty rack for 20 to 300µL: **11983456

***Use empty rack for 100 to 1000µL: **11993456**



Volume Range 0.1 to 5 μ L
11917734



Volume Range 2 to 20 μ L
11943466



Volume Range 30 to 300 μ L
11987714



Volume Range 0.1 to 10 μ L
11903466



Volume Range 2 to 50 μ L
16430875



Volume Range 100 to 1000 μ L
11973466



Volume Range 0.1 to 10 μ L
11913466



Volume Range 10 to 100 μ L
11953466



Volume Range 1 to 5 mL
11927754



Volume Range 0.1 to 20 μ L
11977714



Volume Range 20 to 200 μ L
11963466

SureOne Sterile Filtered Pipette Tips Continued

Volume	Style	Graduated	Tip Length	Color	Low Retention	Packaging	Part code	Fit
20 to 200 μ L	Micropoint	At 10, 50, 100, and 200 μ L	58 mm	Clear	No	Racked	11963466	Universal
					Yes	Racked	16440875	Universal
					No	Reload **	11957724	Universal
30 to 300 μ L	Micropoint	At 10, 50, 100, and 200 μ L	58 mm	Clear	No	Racked	11987714	Universal
					Yes	Racked	16450875	Universal
					No	Reload **	11967724	Universal
100 to 1000 μ L	Micropoint	At 100, 200, 500, and 1000 μ L	88 mm	Clear	No	Racked	11973466	Universal
					Yes	Racked	16400875	Universal
					No	Reload ***	11977724	Universal
1 to 5mL		No	117 mm	Clear	No	Bulk	11927754	Gilson PIPETMAN pipettes



Volume Range 0.1 to 5µL
11997724



Volume Range 0.1 to 10µL
10527014



Volume Range 1 to 200µL

Clear
10678325


Yellow
11983426



Volume Range 0.1 to 10µL
10053014

SureOne Standard Pipette Tips

Volume	Style	Graduated	Tip Length	Color	Low Retention	Packaging	Part code	Fit
0.1 to 5µL	Micropoint Nano	At 2.5µL	31 mm	Clear	No	Bulk	11987724	Universal
					No	Racked	11997724	Universal
					No	Racked Sterile	11907734	Universal
0.1 to 10µL	Micropoint Micro	At 2µL	31.3 mm	Clear	No	Bulk	11933416	Universal
					Yes	Bulk	16440885	Universal
					No	Racked	11953416	Universal
					Yes	Racked	16450885	Universal
					No	Racked Sterile	10053014	Universal
					Yes	Racked Sterile	16470885	Universal
					No	Reload	11973416	Universal
Yes	Reload	16460885	Universal					
0.1 to 10µL	Micropoint Extended Length	At 2, 5, 10, and 20µL	46 mm	Clear	No	Bulk	11983416	Universal
					Yes	Bulk	16400885	Universal
					No	Racked	11588402	Universal
					Yes	Racked	16410885	Universal
					No	Racked Sterile	10527014	Universal
					Yes	Racked Sterile	16430885	Universal
					No	Reload	11967714	Universal
Yes	Reload	16420885	Universal					
1 to 200µL	Beveled	No	51 mm	Clear	No	Bulk	11933426	Universal
					Yes	Bulk	16460875	Universal
					No	Racked	10678325	Universal
					Yes	Racked	16470875	Universal
					No	Racked Sterile	11963426	Universal
					Yes	Racked Sterile	16490875	Universal
					No	Reload	11578412	Universal
					Yes	Reload	16480875	Universal
					No	Bulk	10124314	Universal
					No	Racked	11983426	Universal
Yellow	No	Racked Sterile	11903436	Universal				
	No	Reload	11913436	Universal				



Volume Range 1 to 200µL

Clear
11933436

Yellow
11973436



Volume Range 1 to 200µL

Clear
11953446

Yellow
11913446



Volume Range 5 to 300µL
10003414

SureOne Standard Pipette Tips Continued

Volume	Style	Graduated	Tip Length	Color	Low Retention	Packaging	Part code	Fit	
1 to 200µL	Micropoint Thin Wall	No	50 mm	Clear	No	Bulk	11538422	Universal	
					No	Racked	11933436	Universal	
					No	Racked Sterile	11943436	Universal	
					No	Reload	11953436	Universal	
					No	Bulk	11963436	Universal	
				Yellow	No	Racked	11973436	Universal	
					No	Racked Sterile	11983436	Universal	
					No	Reload	11993436	Universal	
					Clear	No	Bulk	11943446	Universal
						No	Racked	11953446	Universal
No	Racked Sterile	11963446	Universal						
1 to 200µL	Beveled Graduated	At 10, 50, 100µL	50 mm	Yellow	No	Reload	12922521	Universal	
					No	Bulk	11903446	Universal	
					No	Racked	11913446	Universal	
				Clear	No	Racked Sterile	11923446	Universal	
					No	Reload	11933446	Universal	
					No	Bulk	11993446	Universal	
					Yes	Bulk	16480885	Universal	
5 to 300µL	Micropoint	At 10, 50, 100, and 200µL	58 mm	Clear	No	Racked	11903456	Universal	
					Yes	Racked	16490885	Universal	
					No	Racked Sterile	10003414	Universal	
					Yes	Racked Sterile	16410895	Universal	
					No	Reload	11538442	Universal	
					Yes	Reload	16400895	Universal	

*Forward technique only



SureOne Standard Pipette Tips Continued

Volume	Style	Graduated	Tip Length	Color	Low Retention	Packaging	Part code	Fit
100 to 1250µL	Micropoint	At 100, 200, 500, and 1000µL	89 mm	Clear	No	Bulk	11548442	Universal
					Yes	Bulk	16420895	
					No	Racked	11568442	
					Yes	Racked	16430895	
					No	Racked Sterile	10164694	
					Yes	Racked Sterile	16450895	
				Blue	No	Reload	11588442	
					Yes	Reload	16440895	
					No	Bulk	10778535	
					No	Racked	10537014	
					No	Racked Sterile	10492725	
					No	Reload	11963456	
250 to 2.5mL	Standard	No	105 mm	Clear	No	Bulk	11987744	Rainin™ EDP2 Pipettes
500 to 2.5mL	Standard	No	115 mm	Clear	No	Bulk	11997744	Eppendorf and Sartorius™ Biohit™ Pipettes
1 to 5mL	Standard	No	150 mm	Clear	No	Bulk	11648138	Eppendorf and Sartorius Biohit Pipettes
1 to 5mL	Standard	No	147 mm	Clear	No	Bulk	17650627	Fisherbrand Elite, Thermo Scientific F1 & F2
1 to 5mL	Standard	No	117 mm	Clear	No	Bulk	11917754	Gilson PIPETMAN Pipettes
1 to 10mL	Standard	No	161 mm	Clear	No	Bulk	11947754	Gilson PIPETMAN Pipettes
1 to 10mL	Standard	No	152 mm	Clear	No	Bulk	11957754	Fisherbrand Elite, Thermo Scientific F1 & F2

SureOne Empty Racks

- Empty racks are used for filter tip reloads and standard non-filter reload
- Inserts not included

Tip Volume	Packaging	Part Code
10 to 20µL	10/Pack, 50/Case	11973456
20 to 300µL	10/Pack, 50/Case	11983456
1250µL	10/Pack, 50/Case	11993456

SureOne Standard Sterility Claims

- E-beam sterile products are certified to be free of RNase, DNase, DNA, endotoxin, bioburden, PCR Inhibitor, and ATP
- Non-sterile products are certified to be free of RNase, DNase, DNA, pyrogen, and endotoxin.



To order a SureOne Standards Certificate, contact Fisher Scientific Customer Service with the Cat. No. and Lot No.

Fisherbrand Specialty Pipette Tips

Specialty Pipette Tips

- Gel-loading tips ideal for use with agarose gels, 2D, and SDS PAGE protein gels
- Extended-length tips can reach the bottom of larger tubes preventing the pipette from touching the tube walls helping to prevent cross contamination
- Low-retention technology paired with the large orifice diameter of genomic tips provide tools for viscous solution pipetting and help prevent cell shearing when working with large molecules or whole cells

Volume	Style	Tip O.D.	Tip Length	Color	Filtered	Packaging	Part code	Fit
1 to 200µL	Gel Loading	0.58mm	69 mm	Clear	No	Bulk	11927734	universal
					No	Racked	11937734	universal
					No	Racked Sterile	11367801	universal
2 to 20µL	Gel Loading	0.58mm	69 mm	Clear	Yes	Racked Sterile	12992571	universal
10 to 100µL	Gel Loading	0.58mm	69 mm	Clear	Yes	Racked Sterile	11967734	universal
1 to 200µL	Extended Length	N/A	91 mm	Clear	No	Bulk	11977734	universal
						Racked	11987734	universal
						Racked Sterile	11997734	universal



Volume Range 1 to 200µL
11937734



Volume Range 20 to 200µL
12942581



Volume Range 2 to 20µL
12992571



Volume Range 1 to 200µL
11927744



Volume Range 100 to 1000µL
11957744



Volume Range 10 to 100µL
11967734



Volume Range 20 to 200µL
11947744



Volume Range 100 to 1000µL
11977744



Volume Range 1 to 200µL
11987734

Fisherbrand Specialty Pipette Tips Continued

Volume	Style	Tip O.D.	Tip Length	Color	Filtered	Packaging	Part code	Fit
20 to 200µL	Extended Length	N/A	91 mm	Clear	Yes	Racked Sterile	12942581	Universal
1 to 200µL	Large Orifice Low Retention	1.5mm	50 mm	Clear	No	Bulk	11917744	Universal
						Racked	11927744	Universal
						Racked Sterile	11937744	Universal
20 to 200µL	Large Orifice Low Retention	1.5mm	50 mm	Clear	Yes	Racked Sterile	11947744	Universal
100 to 1000µL	Large Orifice, Low Retention Graduated at 100, 200, 500, 1000µL	1.5mm	89 mm	Clear	No	Bulk	12992581	Universal
						Racked	11957744	Universal
						Racked Sterile	11967744	Universal
						Racked Sterile	11977744	Universal

Fisherbrand Filter Tip Reload System

The Fisherbrand Filter Tip Reload System offers the same benefits as other SureOne Pipette Tips. It is the only stackable filter reload system and provides an environmental solution without compromising sterility. Reduce waste and save benchtop and storage space.

1. Offset wells allow tips to nest and stack.
2. Rigid inserts with multichannel format offer a robust non-pliable design to help ensure a secure tip seal.
3. Lot numbered inserts for traceability and easy-to-read filter tip volume and Cat. No.
4. Sterile packaging with outer sleeve to protect tips from contamination during refill.
5. Blister-sleeve guide rails align tip inserts for smooth transfer to the rack.



SureOne Rack Features

1. Flush fit inserts sit flush with rack, allowing easy multichannel pipetting.
2. Recycle logo molded into all rack components allows recycling companies to properly recycle the plastic.
3. Notched clasp designed to prop lid open for autoclave venting.
4. Sure stack lid rim prevents racks from sliding around when stacked and carted around a lab.
5. Safe lock latch allows for one-handed opening and closing ability and won't pop open when dropped.
6. Sure snap design securely holds insert in place, preventing it from being lifted out while pipetting.
7. Flip-lid can be either a hinged or lift-off lid; hinged format does not prevent use of a multichannel pipette.

SureOne Non-Filter Tip Reload Packaging

- Product information displayed on box.
- Dust flap prevents debris from contaminating the product.
- Finger knockouts allow easy access to slide out product.
- Canyon cut package allows for easy tip insert removal.

SureOne Pipette Tips Packaging Options



	Bulk	Raked	Raked Sterile	Reload	Filter Tip Reload	Empty Racks
Packs	1,000/Bag	96 Tips/Tray 10 Trays/Pack	96 Tips/Tray 10 Trays/Pack	96 Tips/Tray 10 Trays/Pack	2 x 480 Tips/Tray* 10 Trays/Pack	10/Pack
Cases	10 Bags/Case	5 Packs/Case	5 Packs/Case	5 Packs/Case	5 Packs/Case	5 Packs/Case

*Filter Tip Reload for 1mL: 56 Tips/Tray; 10 Trays/Pack; 5 Packs/Case

SureOne Filter Tip Reloading Steps

Step 1: Firmly holding the rack below the lowest insert, lift the sleeve up and off the rack, moving the entire sleeve to an empty rack. Align the sleeve knockouts with the latch and rear lid posts.

Step 2: Apply pressure on the middle of the top insert to move the lowermost insert into the empty rack.

NOTE: Always place the plastic sleeve on a rack. The sleeve should not rest directly on the countertop.

Fisherbrand Pipette Tip Standardization

The soft pliable collar of these virgin polypropylene tips can be used with many different pipettes, letting you standardize with fewer products.

- Universal fit
- Certified free of RNase, DNase, DNA, and pyrogens
- Racks nest together, with each serving as the lid for the rack beneath
- Packs include five racks of tips, a base, and a dust cover
- Other sizes and packaging options available



Volume	Tip Length	Color	Filtered	Packaging	Quantity	Part code
0 to 200µL	51 mm	Natural	No	Stacked Racks	1,000/Pack	17725264
		Yellow				17735264